

Original Salter-Style® Cannula

For oxygen flows up to 6 LPM with the Salter-Style® curved facepiece that patients have preferred for over 25 years.



Micro Cannula

With a smaller facepiece for improved comfort and appearance for patients who only require flows up to 3 LPM.



High Flow Cannula

For oxygen flows up to 15 LPM where patients prefer comfort and convenience instead of a mask (includes bright green O₂ supply line).



Smaller Salter-Style® Cannulas

For oxygen flows up to 3 LPM and 6 LPM, and one-piece anatomical design for secure positioning as well as patient comfort and freedom.



16SOFT Cannula

Designed to deliver comfort with soft material in both the tubing and facepiece for a soft-to-the-touch feel throughout.



The softness is everywhere. Curved, comfortable prongs fit gently into patients' nostrils without pressing too hard upon the upper lip. Soft, pliable, stretchable tubing curves around the cheeks and ears, lessening any pressure. And the Salter 16SOFT features a totally redesigned slide bolo to resist twisting and turning.