

## PEDIATRIC OXYGEN MASKS

These Salter masks feature a soft, anatomical design that younger patients will find comfortable, especially for long-term use.



## ADULT OXYGEN MASKS

Salter masks offer 3 methods of attachment for a soft, anatomically shaped form to encourage patient use and acceptance. You'll see a difference.



## PERCENT-O<sub>2</sub>-LOCK® VENTURI OXYGEN SYSTEM

These masks offer an easily adjustable, variable Venturi system delivering accurate concentrations of oxygen (24%, 26%, 28%, 31%, 35%, 40%, and 50%) from a single unit. Concentration levels are easily dialed and locked in place.

