



TRACOE button plus

Silicone decannulation button made of silicone

REF 607-45

Length: ♦ 45 mm
 ID: ♦ 8 mm
 OD: ♦ 11 mm
 Including: ♦ neck strap
 ♦ in sterile package

REF 607-70

Length: ♦ 70 mm
 ID: ♦ 8 mm
 OD: ♦ 11 mm
 Including: ♦ neck strap
 ♦ in sterile package

For the therapy of tracheostomised patients after long-term ventilation – safe airway after decannulation. The sealing is very airtight and dislocation save. The patient is able to cough up the stimulus secretion very fast. With the adjustable neck flange the TRACOE button plus can be positioned optimal to the soft neck tissue.

Caution, not suitable for:

- ♦ laryngectomised patients
- ♦ patients with unstable breathing
- ♦ patients, who primarily breathe through the tracheostoma