

P3

Anaesthesia &
Airway Management



www.p3medical.com

About us

P3 Medical Limited is a British company specialising in the design, manufacture and marketing of a range of high quality medical devices primarily for use within the operating theatre.

Our mission is to provide cost-effective solutions to enhance hospital efficiency, patient care and overall safety in operating theatres in the UK and abroad.

Within our two UK manufacturing sites in Bristol (Headquarters) and Preston, our facilities include 4 modern cleanrooms. We are proud to supply both healthcare organisations and international distributors with British products, manufactured to the highest standards. We currently export to more than 50 countries worldwide.

We are also an established UK distributor for a number of well-known international brands. These products complement our own product ranges and enable us to offer a diverse array of high quality products throughout the healthcare industry.

P3 Medical, established in 2001, is a privately owned company with a wealth of experience in providing excellent service both directly to hospitals and via UK and international supply chains.

Vented Tracheal Tube Introducer

Vented Tracheal Tube Introducer

LATEX FREE

- With lumen and twin distal eyes to allow ventilation
- Single-use, reduces the risk of cross infection
- Supplied with a 15mm connector for quick oxygenation
- Extra tip graduations allow accurate judgement of introducer depth within the trachea
- Highly visible outer tube for ease of identification in an emergency, with string to attach to a trolley
- High-quality, robust outer tube prevents damage in transit and storage
- Individually packed, 10 units per box



| Product code | Description | Size |
|--------------|---------------------------------|--------------|
| V14750 | Vented tracheal tube introducer | 14ch x 750mm |