

Drink-Rite® - Providing **safe, controlled** flow of fluids

Drink-Rite®

The Safer Hydration Delivery System

Key Features

- Available in cup or beaker options each provided with 5cc and 10cc inserts. (Approx. cup holding volume 225cc)
- Suitable for users who cannot tolerate thickened fluids
- Suitable for stage 1 thickened fluids (tested with 'Fresubin thickened Stage 1')
- A safe option for users with dementia or learning difficulties to prevent 'gulping' (drinking too quickly)
- Bite-safe and shatterproof
- Dishwasher safe
- Ideal for domestic as well as nursing home use
- CE marked
- Two flow control inserts with each cup. Collects only 5cc or 10cc of fluid with each tip of the device, to deliver a small amount into the mouth.
- Smooth contoured spout to encourage good lip seal.
- Indented lid design to accommodate the nose. Prevents need to tilt the head, for a safer swallowing position.
- Clear vestibule for visibility of contents.

The Drink-Rite® System provides a controlled limited flow of fluid of 5cc or 10cc into the mouth without the user having to tilt their head.

This is important particularly after stroke or in the care of the elderly. Children with development issues will also benefit from the design features of Drink-Rite.

Drink-Rite®

Ordering information:

OSFSA5301	Comprising of: Cup, 5cc & 10cc insert
OSFSA5302	Comprising of: Beaker (no handles), 5cc & 10cc insert



Drink-Rite® Beaker



The Drink-Rite® Cup and Beaker are supplied with 5cc and 10cc inserts



The Drink-Rite® Cup features wide handles to allow easy hold and grip

Weighted to help **reduce spills** for patients with tremor.

Dysphagia Cup

Weighted Drinking Cup

Key Features

- Suitable for warm and cold fluids and thickened fluids
- Weighted to help reduce spills for patients with tremor
- One piece construction with wide handle to allow extended palm – so gripping not required. The cup is supplied with a detachable base for extra stability
- Hand washing recommended
- Made of durable non toxic polycarbonate material, cup base polyethylene
- CE marked

Kapitex Dysphagia Cup is oval shaped to accommodate nose clearance – which prevents the need to tilt the head whilst drinking.



Dysphagia Cup shown with detachable base for extra stability



Dysphagia Cup

Ordering information:

Supplied with detachable base for extra stability

OSFSA5002	200 ml	Almond
OSFSA5003	200 ml	Clear Green

Available in **3 sizes** in packs of 5

Kapi-Cup®

Drinking Cups

Key Features

- Made of polypropylene with smooth moulded edges, bite-safe and shatterproof
- Suitable for warm and cold fluids and thickened fluids
- Small size (30ml) ideal for use with children
- The medium and large sized cups are supplied with a double handle holder
- A choice of single or double handle holders are available in packs of 2
- Dishwasher safe
- CE marked



Kapi-Cup® shown with blue 80 ml cup

Kapi-Cups have a cut out on the non-drinking side so that they can be tilted without interference by the nose.

This allows the drinker to avoid tilting the head back, thus minimising the chance of liquid entering the respiratory tubes and causing choking.



Kapi-Cup®

Ordering information:

OSFSA5501	Small	Green	30 ml	Pack of 5
OSFSA5502	Medium	Blue	80 ml	Pack of 5
OSFSA5503	Large	White	160 ml	Pack of 5
OSFSA5504	Holder Single Handle	Blue		Pack of 2
OSFSA5505	Holder Double Handle	Blue		Pack of 2

Easy grip mugs featuring various lids with **spouts for different fluid flow**.

Sure Grip Mugs

Easy Grip Drinking Mugs

Key Features

- Lids available with small aperture, large aperture (for thickened fluids) and Thermo lid (heat sensitive lid which changes colour if liquid too hot)
- Flow of liquid can be controlled by sucking on spout or depressing centre of lid
- Suitable for warm and cold fluids
- Each lid features a valve which prevents leakage if the mug is tipped over
- Choice of bright colours – to attract the attention of those with dementia
- Optional straw holder attachment which holds the straw in place on the mug rim
- Dishwasher, microwave and freezer safe
- CE marked



Sure Grip Mugs

Ordering information:

OSFSA5201	200 ml	Mug	Clear
OSFSA5202	200 ml	Mug	Blue
OSFSA5203	200 ml	Mug	Yellow
OSFSA5204	200 ml	Mug	Red
OSFSA5205	Lid - Small Hole and Valve		
OSFSA5206	Lid - Large Hole and Valve		
OSFSA5207	Thermo Lid - Small Hole and Valve		
OSFSA5208	Straw Clips		



Kapitex Healthcare Ltd,
Unit 1, Erivan Park
Sandbeck Way, Wetherby
West Yorkshire LS22 7DN

Tel: +44 (0)1937 580211
Email: sales@kapitex.com

www.kapitex.com

