


TuBo-SA

One-piece Endotracheal Tube Holder

TuBo-SA is a one-piece, latex-free tube holder for securing oral endotracheal tubes, double lumen tubes, laryngeal tubes and laryngeal masks. It is suitable for use in adolescents and adults.



Product characteristics

- One-piece tube holder with soft padding
- Appropriate for all types and sizes of endotracheal tubes

Material

Padding: Soft foam with skin-friendly coating

Taffeta strap: polyester

General: Latex-free, suitable for MRI, non-adhesive

Contact us for sample requests or if you have any questions:



Advantages for users

- Easy and quick to use

Advantages for patients

- Skin-friendly and soft materials
- Protection against decubitus ulcer

REF	Product	Details	Sales Unit
20-20	TuBo-SA	Length 26 cm / for adolescents and adults	60 pieces