

Hand-Fix

Hand and Foot Fixation

Hand-Fix is a softly padded one-piece hand and foot fixation. It is applied as self-protection of patients by preventing unintended manipulation of accesses, drainages or life-sustaining technical devices. Hand-Fix is for patients of all ages.

Product characteristics

- Cuff with taffeta straps
- Available in 5 sizes and in different versions



Material

Cuff: Soft foam with a skin-friendly lamination

Taffeta straps: Polyester

Eyelet: Plastic

Fixing ring: Metal

General: Latex-free, non-adhesive,
suitable for MRI (versions without fixing ring)

Contact us for sample requests or if you have any questions:



Advantages for users

- Easy and quick handling
- Secure hold because of the strong hook and loop fastener and tear-resistant straps

Advantages for patients

- Skin-friendly materials
- Soft and thick padding to protect against decubitus ulcer and injuries

Hand-Fix

Hand and Foot Fixation



Hand-Fix 12 + Hand-Fix 12.1



Hand-Fix 22-30



Hand-Fix 17 + Hand-Fix 17.1



Hand-Fix 25.2

REF	Product	Details	Sales Unit
60-01	Hand-Fix 12	cuff with sewn on straps / width 3 cm, length 12 cm, strap length 65 cm / for neonates, infants and small children	30 pieces
60-01-10	Hand-Fix 12.1	cuff with sewn on straps / width 3 cm, length 12 cm, strap length 75 cm / for neonates, infants and small children	30 pieces
60-08	Hand-Fix 17	cuff with sewn on straps and eyelet / width 4 cm, length 17 cm, strap length 75 cm / for small children and children	30 pices
60-08-10	Hand-Fix 17.1	cuff with sewn on straps and eyelet / width 4 cm, length 17 cm, strap length 85 cm / for small children and children	30 pieces
60-50	Hand-Fix Mix	14 x Hand-Fix 12, 14 x Hand-Fix 17 / for neonates, infants, small children and children	28 pieces
60-61	Hand-Fix 22	cuff with fixing ring and straps / width 7 cm, length 22 cm, strap length 150 cm / for adolescents and adults	16 pieces
60-63	Hand-Fix 25	cuff with fixing ring and straps / width 7 cm, length 25 cm, strap length 150 cm / for adolescents and adults	16 pieces
60-63-01	Hand-Fix 25.1	cuff with fixing ring and sewn on straps / width 7 cm, length 25 cm, strap length 200 cm / for adolescents and adults	16 pieces
60-63-03	Hand-Fix 25.2	cuff with fixing ring and straps / width 7 cm, length 25 cm, strap length 200 cm / for adolescents and adults	16 pieces
60-65	Hand-Fix 30	cuff with fixing ring and straps / width 7 cm, length 25 cm, strap length 150 cm / for adolescents and adults	16 pieces
60-65-01	Hand-Fix 30.2	cuff with fixing ring and straps / width 7 cm, length 30 cm, strap length 200 cm / for adolescents and adults	16 pieces