

Splint-Fix

One-piece Fixing and Positioning Splint

The softly padded Splint-Fix splint is used to immobilize fingers, wrists, ankles and elbows. In case of arterial accesses and peripheral infusions, it ensures a smooth, painless and unhindered flow of the infusions. In the operating theatre, the splint assures optimal performance of blood pressure measurement and blood sampling. Splint-Fix is suitable for patients of all ages.



Product characteristics

- Padded splint with flexible aluminum core
- Sewn on straps with hook and loop fastener

Material

Splint: Soft foam with a skin-friendly cover and aluminum core

Straps: Soft elastic

General: Latex-free, MR conditional, non-adhesive

Contact us for sample requests or if you have any questions:

Advantages for users

- Easy and quick to use
- Easily adjustable splint
- Different sizes allow an individual use

Advantages for patients

- Skin-friendly and soft materials
- Minimization of decubitus ulcer



REF	Product	Details	Sales Unit
60-21	Splint-Fix 9	tapered splint: 9 x 3-4 cm (length x width) / 2 straps: 14 x 1,6 cm (length x width) / for neonates, infants and small children	20 pieces
60-22	Splint-Fix 13	splint: 13 x 4 cm (length x width) / 2 straps: 14, 16 x 2,5 cm (length x width) / for neonates, infants and small children	20 pieces
60-23	Splint-Fix 16	splint: 16 x 5 cm (length x width) / 2 straps: 17, 19 x 2,5 cm (length x width) / for small children and children	20 pieces
60-11	Splint-Fix 22	splint: 22 x 9 cm (length x width) / 2 straps: 30, 40, x 5 cm (length x width) / for children, adolescents and adults	14 pieces
60-13	Splint-Fix 33	splint: 33 x 8,5 cm (length x width), 3 straps: 25, 32 and 42 x 5 cm (length x width) / for children, adolescents and adults	10 pieces