

Sizing Configurations, Kits and Accessories

air-Q is available in a wide range of sizes and kits to meet your needs for hospital and EMS use.



air-Q Disposable Intubating Laryngeal Airways (10/Box)						
	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size	Patient Weight
	1.0	Blue	10-3010	air-Q 1.0 Disposable	4.5	<7 kgs
	1.5	Green	10-3015	air-Q 1.5 Disposable	5.0	7-17 kgs
	2.0	Orange	10-3020	air-Q 2.0 Disposable	5.5	17-30 kgs
	2.5	Yellow	10-3025	air-Q 2.5 Disposable	6.5	30-50 kgs
	3.5	Red	10-3035	air-Q 3.5 Disposable	7.5	50-70 kgs
	4.5	Purple	10-3045	air-Q 4.5 Disposable	8.5	70-100 kgs
air-Q Reusable Intubating Laryngeal Airways (1/Pkg)						
	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size	Patient Weight
	0.5	Pink	10-1001	air-Q 0.5 Reusable	3.5	<4 kgs
	1.0	Blue	10-1010	air-Q 1.0 Reusable	4.5	4-7 kgs
	1.5	Green	10-1015	air-Q 1.5 Reusable	5.0	7-17 kgs
	2.0	Orange	10-1020	air-Q 2.0 Reusable	5.5	17-30 kgs
	2.5	Yellow	10-1025	air-Q 2.5 Reusable	6.5	30-50 kgs
	3.5	Red	10-1035	air-Q 3.5 Reusable	7.5	50-70 kgs
	4.5	Purple	10-1045	air-Q 4.5 Reusable	8.5	70-100 kgs

Sizing Configurations, Kits and Accessories



air-Q Disposable Airways/Lubricant/Syringe (10/Box)

	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size	Patient Weight
	1.5	Green	10-3016	1.5 air-Q Disposable w/Lubricant & 10cc Syringe	4.5	7-17 kgs
	2.0	Orange	10-3021	2.0 air-Q Disposable w/Lubricant & 10cc Syringe	5.5	17-30 kgs
	2.5	Yellow	10-3026	2.5 air-Q Disposable w/Lubricant & 10cc Syringe	6.5	30-50 kgs
	3.5	Red	10-3036	3.5 air-Q Disposable w/Lubricant & 10cc Syringe	7.5	50-70 kgs
	4.5	Purple	10-3046	4.5 air-Q Disposable w/Lubricant & 10cc Syringe	8.5	70-100 kgs

air-Qsp Disposable Intubating Laryngeal Airways (10/Box)

	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size	Patient Weight
	1.0	Blue	10-4010	1.0 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway	4.5	4-7 kgs
	1.5	Green	10-4015	1.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway	5.0	7-17 kgs
	2.0	Orange	10-4020	2.0 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway	5.5	17-30 kgs
	2.5	Yellow	10-4025	2.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway	6.5	30-50 kgs
	3.5	Red	10-4035	3.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway	7.5	50-70 kgs
	4.5	Purple	10-4045	4.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway	8.5	70-100 kgs

Sizing Configurations, Kits and Accessories




air-Qsp Reusable Intubating Laryngeal Airways (1/Pkg)

	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size	Patient Weight
	0.5	Pink	10-6005	0.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	4.0	<4 kgs
	1.0	Blue	10-6010	1.0 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	4.5	4-7 kgs
	1.5	Green	10-6015	1.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	5.0	7-17 kgs
	2.0	Orange	10-6020	2.0 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	5.5	17-30 kgs
	2.5	Yellow	10-6025	2.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	6.5	30-50 kgs
	3.5	Red	10-6035	3.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	7.5	50-70 kgs
	4.5	Purple	10-6045	4.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway	8.5	70-100 kgs

air-Q Blocker Intubating Laryngeal Airways Disposable (10/Box)

	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size
	2.5	Yellow	10-5025	2.5 air-Q Blocker Disposable Airway	6.5
	3.5	Red	10-5035	3.5 air-Q Blocker Disposable Airway	7.5
	4.5	Purple	10-5045	4.5 air-Q Blocker Disposable Airway	8.5

Sizing Configurations, Kits and Accessories




air-Q Blocker Disposable Airways/Lubricant/Syringe (10/box)

	Size	Color	Item	Size/Description	Maximum O.E.T.T. Size
	2.5	Yellow	10-5026	2.5 air-Q Blocker Disposable Airway Kit w/Lubricant & 10cc Syringe with Swivel Elbow.	6.5
	3.5	Red	10-5036	3.5 air-Q Blocker Disposable Airway Kit w/Lubricant & 10cc Syringe with Swivel Elbow.	7.5
	4.5	Purple	10-5046	4.5 air-Q Blocker Disposable Airway Kit w/Lubricant & 10cc Syringe with Swivel Elbow.	8.5

air-Q Blocker Tube (10/box)

	Size	Color	Item	Size/Description	Use With
	2.5	Yellow	10-5027	2.5 air-Q Blocker Tube	2.5 air-Q Blocker Disposable Airway - 10-5025
	3.5	Red	10-5037	3.5 air-Q Blocker Tube	3.5 air-Q Blocker Disposable Airway - 10-5035
	4.5	Purple	10-5047	4.5 air-Q Blocker Tube	4.5 air-Q Blocker Disposable Airway - 10-5045



air-Q® Stylet Guide (10/box)

air-Q Size	Connector Color	Item	Max O.E.T.T.	Patient Weight	Stylet Size
1.0	Blue	10-1006	4.5	4-7 kgs	00 (Small)
1.5	Green	10-1006	5.0	7-17 kgs	00 (Small)
2.0	Orange	10-1005	5.5	17-30 kgs	0 (Medium)
2.5	Yellow	10-1005	6.5	30-50 kgs	0 (Medium)
3.5	Red	10-1004	7.5	50-70 kgs	1 (Large)
4.5	Purple	10-1004	8.5	70-100 kgs	1 (Large)