



Special Protection

Tracheal compresses can be used, for example, to provide padding between the tube and the stoma; they also absorb secretions coming from the stoma. Other products, such as special cleaning sets and brushes, help to maintain the extremely important high standards of hygiene. A well-cared tracheostoma ensures that the patient feels comfortable and reduces the complications.

Bibs, scarves, and roll-necks can help to replace the missing protective and filtration functions of the nose. They protect the airway from dryness, cold and dust but also act as coverings to improve the patient's appearance.



REF 903-A

TRACOE *care* Neck Strap

beige, elastic, consists of one part, individually adjustable, with hook

- Packaging unit: ♦ 1 unit
- Dimensions: ♦ max. length approx. 44 cm
♦ width of neck strap approx. 1.5 cm



REF 903-B

TRACOE *care* Neck Strap

padded, adjustable in length, latex-free, with hook-and-loop fastener

- Packaging unit: ♦ 1 unit
♦ 10 units under **REF 903-B-P**
- Dimensions: ♦ max. length approx. 45 cm
♦ width of neck strap approx. 3.0 cm



REF 903-D

TRACOE *care* Neck Strap

white, wide, padded, elastic, consists of two parts, individually adjustable, with hook-and-loop fastener

- Packaging unit: ♦ 1 unit
- Dimensions: ♦ max. length approx. 49 cm
♦ width of neck strap approx. 2.5 cm


REF 903-E
TRACOE *care* Neck Strap

white, wide, consists of one part, with hook-and-loop fastener

- Packaging unit:
- ◆ 1 unit
 - ◆ 10 units under **REF 903-E-P**
- Dimensions:
- ◆ max. length approx. 42 cm
 - ◆ width of neck strap approx. 3.5 cm


REF 903-F
TRACOE *care* Neck Strap

white, wide, consists of two parts, individually adjustable, with hook-and-loop fastener

- Packaging unit:
- ◆ 1 unit
 - ◆ 10 units under **REF 903-F-P**
- Dimensions:
- ◆ max. length approx. 43 cm
 - ◆ width of neck strap approx. 2.6 cm


REF 903-G

TRACOE *care* Neck Strap

white, padded, consists of two parts, individually adjustable, with hook-and-loop fastener, for neonates, infants and children

- Packaging unit:
- ◆ 1 unit
- Dimensions:
- ◆ max. length approx. 29 cm
 - ◆ width of neck strap approx. 2.6 cm

TRACOE *care* Neck Straps

Ordering Number	max. Length (approx.) cm	Width (approx.) cm	Packaging Unit
REF 903-A	44	1.5	1 unit
REF 903-B	45	3.0	1 unit
REF 903-B-P	45	3.0	10 units
REF 903-D	49	2.5	1 unit
REF 903-E	42	3.5	1 unit
REF 903-E-P	42	3.5	10 units
REF 903-F	43	2.6	1 unit
REF 903-F-P	43	2.6	10 units
REF 903-G	29	2.6	1 unit



TRACOE tube clean Cleaning Kit

with cleaning powder REF 932

REF 930-A

- ◆ with 50 small cleaning swabs (REF 935)

REF 930-B

- ◆ with 50 large cleaning swabs (REF 936)

Including: ◆ tube clean cleaning tub with mesh basket (REF 934)
◆ 1 cleaning brush (REF 938-A)

For cleaning and maintaining tracheostomy tubes, stoma and grid buttons and the corresponding spare parts and accessories.



REF 932

TRACOE tube clean Cleaning Powder

Contents: ◆ 400 g

Cleaning powder for cleaning and maintaining tracheostomy tubes, stoma and grid buttons and the corresponding spare parts and accessories.
approx. 80 applications, approx. 5 g/500ml water



REF 935

TRACOE tube clean Cleaning Swabs small

Packaging unit: ◆ 50 units

Bendable plastic swabs with polyurethane sponge tip for cleaning tracheostomy tubes with an inner diameter of 3 – 7 mm. For single use.



REF 934

TRACOE tube clean Cleaning Tub

with mesh basket

Volume: ◆ 1000ml

For use with the solution, to clean tracheostomy tubes, stoma and grid buttons and the corresponding spare parts and accessories.



TRACOE tube clean Cleaning Kit with concentrated cleaning liquid REF 933

REF 931-A

- ◆ with 50 small cleaning swabs (REF 935)

REF 931-B

- ◆ with 50 large cleaning swabs (REF 936)

Including: ◆ tube clean cleaning tub (REF 934)
◆ 1 cleaning brush (REF 938-A)

For quick cleaning and maintaining tracheostomy tubes, stoma and grid buttons and the corresponding spare parts and accessories.



REF 933

TRACOE tube clean Concentrated Cleaning Liquid

Contents: ◆ 250ml

Liquid cleaning concentrate for cleaning and maintaining tracheostomy tubes, stoma and grid buttons and the corresponding spare parts and accessories. The containing enzymes penetrate and dissolve all organic residues. Rapidly acting.
approx. 100 applications, approx. 2.5 ml/500ml water

REF 936

TRACOE tube clean Cleaning Swabs large

Packaging unit: ◆ 50 units

Bendable plastic swabs with polyurethane sponge tip for cleaning tracheostomy tubes with an inner diameter of 8 – 14mm. For single use.



More information about our cleaning products in our product video on our website
www.tracoe.com/en/media-library



REF 938

TRACOE tube clean Cleaning Brushes with microfibre tip

Packaging unit: ♦ 10 units

Ordering Number	Bristle-diameter mm	for tube Inside Diameter mm	approx. Total Length mm
938-A	7	5/6	210
938-B	9	7/8	210
938-C	11	9/10	210

Angled brush with extra-soft nylon bristles for removing incrustations. To reduce scratching, the brushes have a microfibre tip.



REF 902

TRACOE care Paraffin

Contents: ♦ 25 ml

Lubricates and protects TRACOE *comfort* tubes. Lets inner cannulas be easily inserted and removed, assures best possible fit and good quality of speech.



REF 677

TRACOE care Lubricating Gel

Packaging unit: ♦ 30 Sachets à 2.7 g
♦ sterile

Before placing tracheostomy tubes and/or endotracheal tubes slightly coat them with the lubricating gel. The perfectly dosed gel is colour- and odorless, greaseless, water-soluble and does not stick.



REF 970

TRACOE care Skin Cloth

Packaging unit: ♦ 50 units

For daily cleansing and care around the stoma: pH-balanced skincare products with a pH of 5.0 – 6.0.